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The screenshot displays the GCMS-QP2010 software interface. On the left is a vertical toolbar with icons for Quantitation, Calibration Curve, Batch Processing, Data Comparison, Report Format, Library Editor, and Guide. The main window is divided into several panels:

- Parameters Table:**

Index	Parameter	Upper/Lower	# of Hit
2	Mol Wt	284-500	No need to set 55960
3	Capd Name	<input checked="" type="checkbox"/> Match Case	55960
- Hit List Table:**


Hit#	Capd Name	Mol Wt	Formula
1	Propane, 1,1,1,2,2,3-hexachloro-3,3-difluoro- δ	284	C9Cl6F2
2	3,4-Dibromothiophene-2-carboxylic acid $\delta\delta$ 3,4-D	284	C5H2Br2O2S
- Mass Spectrum:** A plot of Intensity (x100) vs. m/z. The base peak is at m/z 85. Other significant peaks are at m/z 117, 121, 129, 166, 179, 201, 231, 251, and 267. A chemical structure of Propane, 1,1,1,2,2,3-hexachloro-3,3-difluoro- δ is shown above the spectrum.
- Search Details:**

CAS# : 661-96-1 Mol Wt 284 Serial# 85322
Capd Name : Propane, 1,1,1,2,2,3-hexachloro-3,3-difluoro- $\delta\delta$ 1,1-Difluoro-1,2,2,3,3,3-hexachloro- δ
Capd Form : C9Cl6F2
Class : No Class Flags.

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The screenshot shows the GCMS Real Time Analysis (Admin) software interface. The main window displays a table with the following data:


	Vial#	Sample Name	Sample ID	Sample Type	Analysis Type	Method File	Data File	Level#	Inj. Volume
1	1		VNK-0001	0: Unknown	IT QT	n\Sample\SIM.qm	t_std_001.qgd	1	1
2	2			0: Unknown	IT QT	n\Sample\SIM.qm	t_std_002.qgd	1	1
3	3			0: Unknown	IT QT	n\Sample\SIM.qm	t_std_003.qgd	1	1
4				0: Unknown	IT QT	n\Sample\SIM.qm	known_001.qgd	1	1
5	1			0: Unknown	IT QT			1	1

On the left side of the interface, there is a vertical toolbar with several icons. A circled number "2" with an arrow points to the "Start" button, which is represented by a green play button icon.

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- ③ 输入 Wizard 表 1/7


The screenshot shows the "Compound Table Wizard 1/7" dialog box. The window title is "Compound Table Wizard 1/7". On the left, there is a graphic of a molecular structure and a table. The main text area contains the following text: "This wizard automatically creates a Compound Table from the spectra registered in the Spectrum Processing Table. Keep in mind that the previous Compound Table is deleted." Below this text, there are two radio button options: "Use current Spectrum Process Table" (which is selected) and "Integration of TIC". Below these options are two buttons: "Integration Param..." and "Spectrum Format...". At the bottom of the dialog, there are four buttons: "< 上一步(B)", "下一步(N) >", "取消", and "帮助".

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
Zgro: Not Forced

Weighted: None

Format of Concentration: Decimal Signific: 5

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Compound Table Wizard 5/7

First, enter the standard concentration for each level. Then set the amount of internal standard to use in the Internal Standard field. In the Number of Reference Ions field, enter zero to use no reference ion.

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