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按下 **Calculate** 按钮

Isocratic Separations

Select your existing HPLC conditions here:

Column Length	Column Diameter	Particle Size	MW of Sample	Flow Rate	Temperature	Retention Factor
15 cm	3.9 mm	4 μm	100	0.60	40	12
12.5 cm	3.0 mm	3.5 μm	200	[mL/min]	[Celsius]	[k]
10 cm	2.1 mm	3 μm	300			

Solvent Composition

%Aqueous: 70 MeOH: 30 MeCN: 30

Pressure Units: bar PSI MPa

Injection Volume: 10 [μL]

Results With Existing Column

N: 9464 Plates

Run Time: 10.2 min

Pressure: 1164 PSI

Select Maximum Column Pressure 10000 PSI **Change**

Click here for your ACQUITY UPLC™ column choices.

Calculate Print Reset Exit

Results | Explanation of Results

Conditions of Equal Efficiency			Conditions of Maximum Efficiency (10 cm Column Length)			Conditions of Shortest Analysis Time (5 cm Column Length)		
Column Length	5	cm	Max. Plates	21160	Run Time	1.4	min	Plates
Run Time	2.5	min	Run Time	11.5	min	Plates	7777	
Flow Rate	0.600	mL/min	Flow Rate	0.264	mL/min	Flow Rate	1.075	mL/min
Pressure	5565	PSI	Pressure	4836	PSI	Pressure	10000	PSI
Column ID	2.1	mm	Column ID	2.1	mm	Column ID	2.1	mm
Inj. Volume	2.5	μL	Inj. Volume	4.9	μL	Inj. Volume	2.5	μL

maximum value = 200. The retention factor is the analysis time divided by the column dead time minus 1.

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Conditions of Equal Efficiency		Conditions of Maximum Efficiency (10 cm Column Length)		Conditions of Shortest Analysis Time (5 cm Column Length)	
Column Length	5 cm	Max. Plates	21160	Run Time	1.4 min
Run Time	2.5 min	Run Time	11.5 min	Plates	7777
Flow Rate	0.600 mL/min	Flow Rate	0.264 mL/min	Flow Rate	1.075 mL/min
Pressure*	5565 PSI	Pressure*	4836 PSI	Pressure*	10000 PSI
Column ID	2.1 mm	Column ID	2.1 mm	Column ID	2.1 mm
Inj. Volume	2.5 µL	Inj. Volume	4.9 µL	Inj. Volume	2.5 µL

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The screenshot shows the 'Gradient Separations' software window. It includes a menu bar (File, Run, Help), a title bar, and a main content area. The main area is divided into several sections: 'Select your existing HPLC conditions here:' with a table of parameters (Column Length, Column Diameter, Particle Size, MW of Sample, Flow Rate, Temperature, Injection Volume, Instrument Delay); 'Organic Modifier in B' (MeOH %100); 'Pressure Units' (bar, PSI, MPa); 'Results With Existing Column' (Peak Capacity, Pressure); and a 'Gradient Table' with columns for Step, Time, Flow, ΔA , ΔB , Time Segment (min), and Column Volumes. A 'Calculate' button is at the bottom right. Callouts point to various elements: '现有的 HPLC 条件' (Existing HPLC conditions) points to the parameter table; '现有 HPLC 的梯度表' (Existing HPLC gradient table) points to the gradient table; '调整 UPLC 色谱柱的最大反压' (Adjust the maximum backpressure of the UPLC column) points to the 'Select Maximum Column Pressure' field; '计算 (估计) 出的 HPLC 结果' (Calculated (estimated) HPLC results) points to the 'Results With Existing Column' section; '计算出的时间段和柱体积' (Calculated time segments and column volumes) points to the 'Time Segment' and 'Column Volumes' columns; and '屏幕帮助窗口' (Screen help window) points to the bottom status bar.

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Peak Capacity:	118	Column Length:	15 cm	Particle Size:	5 μm
Pressure (PSI):	1455	Column Diameter:	4.6 mm	MW of Sample:	200

Length [cm]	ID [mm]	Flow mL/min	Peak Capacity	Run Time [min]	Pressure [PSI]	Injection Volume [μL]	Detailed Gradient Profile
5 cm	2.1	0.919	---	---	---	---	View
10 cm	2.1	0.919	---	---	---	---	View

Length [cm]	ID [mm]	Flow mL/min	Peak Capacity	Run Time [min]	Pressure [PSI]	Injection Volume [μL]	Detailed Gradient Profile
5 cm	1.0	0.071	65	5.1	3768	0.2	View
10 cm	1.0	0.071	116	10.1	7536	0.3	View

--- indicates maximum UPLC™ pressure limit exceeded. Please increase UPLC™ maximum pressure, lower HPLC flow rate and/or increase temperature.

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Length [cm]	ID [mm]	Flow mL/min	Peak Capacity	Run Time [min]	Pressure [PSI]	Injection Volume [μL]	Detailed Gradient Profile
5 cm	2.1	0.574	179	15.0	6319	0.7	View
10 cm	2.1	0.510	227	15.0	12501	1.4	View

Length [cm]	ID [mm]	Flow mL/min	Peak Capacity	Run Time [min]	Pressure [PSI]	Injection Volume [μL]	Detailed Gradient Profile
5 cm	2.1	0.998	119	2.5	11968	0.7	View
10 cm	2.1	0.532	119	2.9	12501	1.4	View

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© 2005 Waters Corporation