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Determination of Benzoate in Liquid Food Products by Reagent-Free Ion Chromatography

Yan Jing

(Shanghai Institute of measuring and testing, Physical & Chemical Analysis Department, Shanghai, 201203, China)

Abstract A simple ion chromatography method was described for the direct determination of benzoate in liquid food products. Brand-new ICS-3000 ion chromatography system with reagent free eluent generation device were used with high capacity IonPac AS18 column by which standard solutions and actual samples were analysed and the results are satisfied. Only dilution and de-gas process need to be done before injection. This method which has the advantages of easy operation, good reproducibility, high sensitivity etc and the spiked results are very good as well. Proven by our cooperation labs, US FDA and US Dionex, this method can be used as the standard method to analyse Benzoate in food and Beverage matrix.

Key words Ion chromatography; food; Benzoate

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- [1] Velsicol Chemical Corporation, Product Information Bulletin; Sodium Benzoate Food Grade, FCC; <http://www.velsicol.com/aspa/documents/pibs/pdfpib42.pdf>
- [2] U. S. Food and Drug Administration. Code of Federal Regulations; Title 21 Volume 6, Section 582.3733 sodium benzoate; U. S. Government Printing Office, Washington, DC, Revised April 1, 2003
- [3] European Parliament and Council Directive No. 95/2/EC. Food Additives Other Than Colours and Sweeteners; February, 1995
- [4] Association of Official Analytical Chemists (AOAC). Benzoic Acid and Sorbic Acid in Food; Method 983.16; Gaithersburg, MD, 1995