



核苷酸的分离

Sunnest RP-AQUA 3 μ m150 x 4.6 mm i.d.色谱柱: Sunnest RP-AQUA, 3 μ m 150 x 4.6 mm

仪器: LC-5090

流动相: 20mM Phosphate buffer pH6.0

流速: 1.0 mL/min

柱温: 40 °C

检测器: UV@250nm

样品: 1 = 5'-GDP, 2 = 5'-GMP, 3 = 5'-ATP, 4 = 5'-ADP, 5 = 5'-AMP

