


氨基酸(LC/MS)的分离

氨基酸(LC/MS)

Sunniest RP-AQUA 5 μ m, 150 x 2.1 mm i.d.

色谱柱: Sunniest RP-AQUA 5 μ m, 150 x 2.0 mm

流动相:

A) 5mM HFBA (Heptafluorobutyric acid)

B) 5mM HFBA in Acetonitrile/water(9/1)

%B 0% to 20% in 20 min

流速: 0.2mL/min

柱温: 40 $^{\circ}$ C

检测器: Quattro Micro API (ESI positive) SIM and UV@210nm