


缬沙坦及有关物质的分离


Sunniest C18 3 μ m 150 x 4.6 mm i.d.

Sunniest C18 3 μ m 150 x 4.6 mm i.d.



1 Unknown

2



3 Unknown

色谱柱: Sunniest C18, 3 μ m 150 x 4.6 mm

流动相: Acetonitrile/Water/Glacial Acetic Acid = 500/500/1

流速: 0.4 mL/min

柱温: 30 $^{\circ}$ C

检测器: UV@225nm

测试探头: Valsartan 0.5mg/mL in mobile phase

进样量: 10 μ L